

## Protect patients and staff by staying on top of healthcare hygiene standards

Trusted hygiene solutions for hospitals



Tork is your trusted partner in maintaining infection prevention and hygiene compliance within your hospital. Tork solutions fulfil strict healthcare requirements and fit into the workflow of health professionals and cleaners. The design of Tork dispensers is functional and attractive, signalling to visitors and staff that hygiene is taken seriously. Overall, our innovative solutions can help you meet hygiene and cleaning standards, all while meeting patient and staff expectations.



### Enhance your hygiene systems and outcomes

A clean and hygienic environment plays a vital role in providing quality healthcare and positive outcomes. You can depend on Tork hygiene solutions in the care environment.

- Tork PeakServe® hand towel dispensers help reduce time restocking dispensers and increase time for hospital staff to focus on maintaining a clean and safe environment
- Convenient placement of hand sanitiser encourages hand hygiene to reduce the spread of infections
- Smooth rounded edges on dispensers are easy to wipe down and clean
- Tork provides additional support with dispenser placement, tips, Tork VR hand hygiene training and hand hygiene signage



### Improve perceptions and the healthcare experience

Patient and visitor perceptions of care are influenced by many factors. Tork solutions help make hygiene a positive factor.

- Compact, whisper-quiet hand towel dispensers are perfect for use in patient rooms and quiet zones
- Increase accessibility with the silent and slim Tork PeakServe hand towel dispensers, featuring the lowest consistent pull force of any high-capacity hand wiping system
- Ensure public and patient care areas are clean and welcoming with functional and attractive Tork dispensers



### Enhance cleaning efficiency and staff productivity

A clean and hygienic environment is a priority in hospitals. Tork solutions help ensure cleaning time is productive and done where it's most needed.

- Tork EasyCube® connected dispensers provide real-time information to cleaners to help them act on what's needed and where – freeing up more time for hygiene-critical tasks with less time checking and refilling dispensers
- Tork PeakServe high-capacity dispensers shift more time to hygiene-critical tasks because they use 50% compressed towels that allow more towels on the trolley and lead to fewer restocking trips
- All Tork dispensers are easy to load, freeing up more time for cleaning
- Low-maintenance systems enhance staff efficiency

# Recommended Tork solutions



## Patient and Treatment Rooms

Tork user-friendly systems make it easy to prevent cross-contamination while improving image and the experience.



552100  
Tork Xpress® Mini Multifold Hand Towel Dispenser



560000  
Tork Liquid and Spray Skincare Dispenser



420110  
Tork Alcohol Liquid Hand Sanitiser



682000  
Tork SmartOne® Twin Mini Toilet Roll Dispenser



473264  
Tork Reflex Wiping Paper Plus Centrefeed Roll



140278  
Tork Extra Soft Facial Tissue Cube



140280  
Tork Extra Soft Facial Tissues



## Public Washrooms

Tork user-friendly systems make it easy to prevent cross-contamination while improving image and experience. Tork PeakServe high-capacity dispensers also help free up more time for hygiene-critical tasks with less time checking and refilling dispensers.



552500  
Tork PeakServe® Continuous™ Hand Towel Dispenser



551000  
Tork Matic® Hand Towel Roll Dispenser



571600  
Tork Foam Skincare Automatic Dispenser



420701  
Tork Extra Mild Foam Soap



682000  
Tork SmartOne Twin Mini Toilet Roll Dispenser



140278  
Tork Extra Soft Facial Tissue Cube



140280  
Tork Extra Soft Facial Tissues



## Common Areas

Tork provides high-quality refills with innovative dispensers to create a safe and welcoming atmosphere to enhance patient safety and visitor satisfaction.



420103  
Tork Alcohol Gel Hand Sanitiser



420101  
Tork Alcohol Foam Hand Sanitiser



511055  
Tork Hygiene Floor Stand



571600  
Tork Foam Skincare Automatic Dispenser




140278  
Tork Extra Soft Facial Tissue Cube



140280  
Tork Extra Soft Facial Tissues

## Recommended Tork solutions

Continued



### Cafeteria

Tork solutions ensure the very best practices for cleanliness and comfort in these busy gathering spots.



272611  
Tork Xpressnap® Tabletop  
Napkin Dispenser



272901  
Tork Xpressnap Fit® Countertop  
Napkin Dispenser



### Cleaning

Tork solutions support the strictest hygiene protocols to help prevent infections and improve safety.



654000  
Tork Folded Wiper/Cloth Dispenser



530178  
Tork Heavy-Duty Cleaning Cloth

How else can we help?  
Get in touch to talk about  
your business.

tork.co.uk or tork.ie  
+44 (0)1582 677570 or  
+353 (0)1 7930 150  
TorkCS.uk@Essity.com or  
TorkCS.ie@Essity.com