



Think ahead.

Tork Skincare

Meet high hygiene expectations, with Tork skincare solutions



We have 50 years of experience

As the global leading professional hygiene brand, Tork has the solutions and innovations companies trust to help them secure the new hygiene standard.

In addition to our innovative portfolio, we have a broad, effective assortment of soap, sanitizer and dispensers.

Not only are Tork skincare solutions effective, the award-winning dispenser designs coordinate aesthetically with the hand towel and toilet paper systems customers are already using.

What's more, Tork skincare dispensers are certified easy-to-use,¹ leaving a positive impression on guests, encouraging hand hygiene compliance and reducing cross-contamination.

Customers already depend on Tork for their paper hand towel and toilet tissue needs and look to Tork for quality and effective skincare solutions that address both new hygiene expectations and ongoing business challenges.

With the Tork skincare assortment, staff can focus on delivering great customer service – not addressing washroom problems.



¹ Certified by the Swedish Rheumatism association, tested according to The Design for All Test, published in the Swedish Design Research Journal, 2011-05-20

High-quality hand hygiene solutions

A comprehensive product selection:

A complete range of effective sanitizers and soaps with carefully chosen ingredients.

Flexible manual or sensor dispensers help secure a higher hygiene standard.

Tork Hygiene Stands keep hand hygiene top of mind and accessible for employees and customers when no water is available.

Securing best-in-class hand hygiene

Meeting high standards for efficacy: With a variety of antimicrobial soap and hand sanitizers in liquid, gel and foam form – all of which are tested to either an EN1499 or EN1500 standard.

Hygiene for all: Tork skincare dispensers' soft push buttons are certified "Easy to Use" by the Swedish Rheumatism Association,³ leading to easier dispensing for children and the elderly.

A hygienic system: Factory-sealed bottles, with a new pump for every refill, help reduce risk of contamination, compared to open bulk systems where 25% of open soap dispensers can be contaminated.⁴

44%

avoid washrooms due to hygiene concerns²



² Based on results from Hygiene Matters survey 2018 of 15 530 respondents from 15 markets

³ Certified by the Swedish Rheumatism association, tested according to The Design for All Test, published in the Swedish Design Research Journal, 2011-05-20

⁴ Bulk Soap Contamination Study, Dr. Charles P. Gerba, University of Arizona, 2007

Maximize handwashing efficiency

Employees need to practice hand hygiene often and effectively to prevent germ spread and reduce cross contamination, so they can maximize their time spent working.

To ensure washroom and handwashing efficiency, cleaners need to be able to refill dispensers quickly and correctly, so they can focus on other pressing cleaning tasks.



Easy-to-use: Tork dispensers are “Easy to Use,”⁶ and take less than 10 seconds to refill,⁷ saving time and energy for cleaning staff. Tork also provides a suite of washroom hygiene products including toilet tissue, hand towels, haircare and skincare. Each dispenser comes with one common Tork key, which makes refilling more efficient.

Reduces consumption: Measured, consistent dispensing reduces consumption by up to 50%,⁸ meaning less refills and increased soap availability.⁹

High capacity: Tork Foam Soap refills can serve up to 1,650 visitors, so less frequent refilling is required.⁹

⁵ Based on the documented results achieved by three Tork EasyCube customers, measured before and after the implementation of Tork EasyCube®

⁶ Certified by the Swedish Rheumatism association, tested according to The Design for All Test, published in the Swedish Design Research Journal, 2011-05-20

⁷ Based on internal panel test soap refilling time

⁸ Essity usage studies on 20,000 hand washes, comparing different liquid soaps (including bulk soap) to Tork Foam and Spray Soap Systems

⁹ TMTC Consumption Test Essity usage studies on 20,000 hand washes, comparing different liquid soaps (including bulk soap) to Tork Foam and Spray Soap Systems

Functional and intuitive, Tork skincare dispensers limit hassle

Robust, hassle-free dispensers: Tork dispensers are proven to perform hassle-free up to 1MM pushes.¹⁰

Skincare, anywhere: Tork skincare dispensers come with a variety of dispensing mechanisms to best suit facility-specific needs so that hand hygiene can happen anywhere.

Hygiene inclusivity for all: Tork skincare dispensers' soft push buttons are certified "Easy to Use"¹¹ by the Swedish Rheumatism Association leading to easier dispensing and encouraging hand hygiene compliance.



Easy-to-maintain:

- Tork dispensers take less than 10 seconds to refill, saving time and energy for cleaning staff.¹²
- Tork provides a suite of washroom hygiene products including toilet tissue, hand towels, haircare and skincare. Each dispenser comes with one common Tork key, which makes refilling more efficient.

"Easy to Use" certified:¹¹ leading to easier dispensing

High-capacity: serving up to 1,650 visitors¹³

Touch free: ensuring measured, consistent dispensing



Factory-sealed bottle: helping reduce risk of cross contamination compared to open bulk systems

Collapsible bottle: leading to 70% less waste volume¹⁰

¹⁰ Based on Essity test.

¹¹ Certified by the Swedish Rheumatism association, tested according to The Design for All Test, published in the Swedish Design Research Journal, 2011

¹² Based on internal panel test soap refilling time

¹³ TMTC Consumption Test Essity usage studies on 20,000 hand washes, comparing different liquid soaps (including bulk soap) to Tork Foam and Spray Soap Systems

Tork skincare, delivering premium hand washing experiences for 50+ years

Washroom aesthetics are an important consideration for decision makers.

By creating a pleasant environment, you can generate a lasting, positive impression for both guests and employees.



Award-winning:

Two-time Red Dot Design award-winning Tork dispensers are recognized for their sleek look and feel.

Designed with the whole washroom in mind, Tork dispensers deliver a premium guest experience, and have coordinated designs (Elevation and Image Design) across all dispenser types.



Skincare soap and lotion complete the experience:

High-end Tork formulas provide visitors an extra touch of luxury.

Tork skincare delivers efficacy and sustainability

Companies value sustainable products and solutions to meet sustainability business targets. In tandem, employees feel more positively about their workplace when sustainable products are implemented.

With Tork skincare formulations, decision makers do not have to compromise on quality to meet their sustainability targets.



Nearly all Tork skincare refills have third-party EU eco-label certification, helping to reduce environmental impact.



¹⁴ Based on Essity test

¹⁵ Proven by third part lab test according to OECD301B

¹⁶ Essity usage studies on 20,000 hand washes, comparing different liquid soaps (including basic soap) to Tork foam and Spray Soap Systems

¹⁷ Based on an internal panel test comparing Tork Mild Foam Soap

Specialized hand hygiene training:

Ella's Handwashing

Adventure teaches children hand hygiene, building critical habits at a young age.



The award-winning Tork VR Clean Hands Training, which is also available on desktop, promotes hand hygiene compliance in acute and long-term care facilities. This new training method is a result of a collaboration between Tork, university hospitals and world-leading hand hygiene experts.



Tork skincare, developed with almost every situation in mind



To meet your sustainability needs



Tork Clarity Foam Soap
SKU: 520201



Tough on viruses and bacteria. Gentle on everything else.



Tork Antimicrobial Liquid Soap and Alcohol-Free Foam Hand Sanitizer
SKU: 420709-424709-420710-424710-520202-520402



For when proper hand sanitation is critical



Tork Alcohol Liquid Hand Sanitiser
SKU: 420114
adapt to local SKU number



For the most demanding environments



Tork Oil and Grease Liquid Soap
SKU: 420401



For professional kitchens



Tork Odour-Control Hand Wash
SKU: 424011



To deliver a premium handwashing experience



Tork Luxury Soft Liquid Soap
SKU: 420901



To clean and soften both hair and skin



Tork Shower Cream
SKU: 420601



For those with sensitive skin



Tork Extra Mild Liquid Soap
SKU: 420701



To combat germs, while remaining gentle on skin



Tork Alcohol Gel Hand Sanitizer
SKU: 424103
adapt to local SKU number

The full Tork skincare offer for an effective handwashing experience

With our full range of formulations and dispensers for all washroom needs, we can guide you to select the right hand hygiene solution for your washroom.

Some SKU numbers must be adapted locally.



Tork Skincare Foam Soap and Liquid Soap

 S4, S3



- Dispensing system easy to use for children and elderly
- Sealed system with collapsible bottle improves hygiene and reduces waste

Dispensers

Art. No	Name	Design Line	Refill ID	Colour	Material	Size HxWxD (mm)	CON/ TRP		
561600	Tork Skincare Dispenser – with Intuition™ sensor	Elevation	S4	White	ABS	278x112x128	6	●	●
561608	Tork Skincare Dispenser – with Intuition™ sensor	Elevation	S4	Black	ABS	278x112x128	6	●	●
561500	Tork Skincare Dispenser	Elevation	S4	White	ABS	292x113x114	12		●
561508	Tork Skincare Dispenser	Elevation	S4	Black	ABS	292x113x114	12		●
460009	Tork Skincare Dispenser - with Intuition™ sensor	Image Design	S4	Metallic/Black	Stainless Steel/ABS	278x116x130	1	●	●
460010	Tork Skincare Dispenser	Image Design	S4	Metallic/Black	Stainless Steel/ABS	289x106x107	1		●




Cosmetic Refills

Art. No	Name	Refill ID	Quality	Fragrance	Colour	Benefit	Content/ Unit (ml)	No. of Doses	CON/ TRP			
520901	Tork Luxury Soft Foam Soap	S4	Premium	Sophisticated Floral Scent	Dye Free	Rich & Creamy Foam Handwash	1 000	2 500	6		●	
New 520201	Tork Clarity Hand Washing Foam Soap	S4	Premium	Perfume Free	Dye Free	99% of ingredients from natural origin	1 000	2 500	6	●	●	
520701	Tork Extra Mild Foam Soap	S4	Premium	Perfume Free	Dye Free	Frequent Handwash	1 000	2 500	6	●	●	●
520501	Tork Mild Foam Soap	S4	Premium	Fresh Scent	Dye Free	Regular Handwash	1 000	2 500	6	●	●	

Tork Luxury Soft Foam Soap S3 fit in foam soap dispensers from Tork Aluminium Line 453000 that are already installed.

Antimicrobial (Biocide) Refills¹

Art. No	Name	Refill ID	Quality	Fragrance	Colour	Efficacy	Efficacy Norms	Content/ Unit (ml)	Dose-Size (ml)	CON/ TRP	
520800	Tork Antimicrobial Foam Soap Improved Formulation	S4	Premium	Perfume free	Dye free	Bacteria Yeast & Fungi Viruses (enveloped)	EN 1499, EN 13727 EN 13624, EN 14476 (all enveloped viruses)	1 000	0.6	6	●

¹ Use biocides safely. Always read the label and product information before use.



424011



424709

Cosmetic Refills

Art. No	Name	Refill ID	Quality	Fragrance	Colour	Benefit	Content/ Unit (ml)	No. of Doses	CON/ TRP		Other Eco-labels		
New 424011	Tork Odor-Control Hand Washing Liquid Soap	S4	Premium	Perfume Free	Dye Free	Odor control, Food Preparation, Handwash	1 000	1 000	6				
New 424709-424710	Tork Antimicrobial Liquid Soap	S4	Premium	Perfume Free	Dye Free	Tough on viruses and bacteria. Gentle on everything else	1 000	1 000	6				

Tork Skincare Mini Liquid Soap

S2


561000



561008

- Compact dispenser size fits even in the smallest washroom
- Dispensing system easy to use for children and elderly
- Sealed system with collapsible bottle, ensures good hand hygiene

Dispensers

Art. No	Name	Design Line	Refill ID	Colour	Material	Size HxWxD (mm)	CON/ TRP	
561000	Tork Mini Liquid Skincare Dispenser	Elevation	S2	White	ABS	206x112x114	12	
561008	Tork Mini Liquid Skincare Dispenser	Elevation	S2	Black	ABS	206x112x114	12	



420502



420602



420202



420702

Cosmetic Refills

Art. No	Name	Refill ID	Quality	Fragrance	Colour	Benefit	Content/ Unit (ml)	No. of Doses	CON/ TRP		Other Eco-labels		
420502	Tork Mild Mini Liquid Soap	S2	Premium	Fresh Scent	Pearl White	Regular Handwash	475	475	8				
420602	Tork Mini Shower Cream	S2	Premium	Fresh Unisex Scent	Blue	Hair & Body	475	475	8				
420202	Tork Non-Perfumed Hand & Body Lotion	S2	Premium	Perfume Free	White	Re-Moisturising Lotion	475	475	8		Swan		
420702	Tork Extra Mild Mini Liquid Soap	S2	Premium	Fresh Unisex Scent	Transparent	Frequent Handwash	475	475	8				

Tork Toilet Seat Cleaner

S2



420302

- For extra hygiene and comfort inside the cubical
- Cleaning for an extra-hygienic toilet seat
- Sealed bottle with single-use pump reduces the risk of cross contamination

Art. No	Name	Refill ID	Quality	Colour	Type	Benefit	Content/Unit (ml)	No. of Doses	CON/TRP	
420302	Tork Toilet Seat Cleaner	S2	Premium	Transparent	Detergent	Extra Hygienic Toilet Seat	475	475	8	

Tork Skincare Liquid Soap

S1/S11



560000

560008

560101

560109

252040

960000

- Provides a large dose of cleaning ingredients
- Sealed system with collapsible bottle improves hygiene and reduces waste
- Dispensing system easy to use for children and elderly

Dispensers

Art. No	Name	Design Line	Refill ID	Colour	Material	Size HxWxD (mm)	CON/TRP	
560000	Tork Liquid and Spray Skincare Dispenser	Elevation	S1/S11	White	ABS	296x112x114	12	
560008	Tork Liquid and Spray Skincare Dispenser	Elevation	S1/S11	Black	ABS	296x112x114	12	
560101	Tork Liquid Skincare Dispenser – Arm Lever	Elevation	S1	White	ABS	296x112x114	8	
560109	Tork Liquid Skincare Dispenser – Arm Lever	Elevation	S1	Black	ABS	296x112x114	8	
252040	Tork Liquid Soap Dispenser		S1	White	Metal	291x113x98	10	
960000	Tork Starter Pack Liquid Soap (Art. No 560000, Art. No 420501)	Elevation	S1	White	ABS	296x112x114	1	



420709

420901

420701

420501

420601

420810

420401

Cosmetic Refills

Art. No	Name	Refill ID	Quality	Fragrance	Colour	Benefit	Content/Unit (ml)	No. of Doses	CON/TRP		Other Eco-labels		
Now 420709-420710	Tork Antimicrobial Liquid Soap	S1	Premium	Perfume Free	Dye Free	Tough on viruses and bacteria. Gentle on everything else	1 000	1 000	6				
420901	Tork Luxury Soft Liquid Soap	S1	Premium	Luxury Jasmine Perfume	Lilac	Rich & Creamy Handwash	1 000	1 000	6		Swan		
420701	Tork Extra Mild Liquid Soap	S1	Premium	Perfume Free	Dye Free	Frequent Handwash	1 000	1 000	6				
420501	Tork Mild Liquid Soap	S1	Premium	Fresh Scent	Pearl White	Regular Handwash	1 000	1 000	6				
420601	Tork Shower Cream	S1	Premium	Fresh Unisex Scent	Blue	Hair & Body	1 000	1 000	6				
420810	Tork Hand Washing Liquid Soap	S1	Premium	Perfume Free	Dye Free	Food Preparation Handwash	1 000	1 000	6				
420401	Tork Oil and Grease Liquid Soap	S1	Premium	Perfume Free	Dye Free	Oil & Grease Handwash	1 000	1 000	6		Swan		



620701

620501

- Functional, high capacity for best in class cost in use
- Reduces Soap usage by up to 50%¹

Skincare Spray Soap Refills

Art. No	Name	Refill ID	Quality	Fragrance	Colour	Benefit	Content/ Unit (ml)	No. of Doses	CON/ TRP		
620701	Tork Spray Soap Non-Perfumed	S11	Universal	Perfume Free	Dye Free	Frequent Handwash	1 000	3 300	6	●	●
620501	Tork Spray Soap	S11	Universal	Fresh Scent	Dye Free	Regular Handwash	1 000	3 300	6	●	●

¹ Essity usage studies on 20.000 hand washes, comparing different liquid soaps (including bulk soap) to Tork Foam and Spray Soap systems.

Tork Alcohol Hand Sanitisers

S1/S2/S4, Freestanding Bottles



420105

420142

424105

511103

590103

914103

Refills - Gel¹

Art.No	Name	Refill ID	Quality	Fragrance	Colour	Efficacy	Efficacy Norms	Content/ Unit (ml)	Dose- Size (ml)	CON/ TRP		
420105	Tork Alcohol Gel Hand Sanitiser	S1	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Mycobacteria Virus	EN 1500, EN 12791, EN 13727, EN 1276, EN 1275, EN 13624, EN 1650, EN 14348, EN 14476 (all enveloped and Rotavirus & Norovirus)	1 000	1	6	●	●
420142	Tork Alcohol Gel Hand Sanitiser	S2	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Mycobacteria Virus	EN 1500, EN 12791, EN 13727, EN 1276, EN 1275, EN 13624, EN 1650, EN 14348, EN 14476 (all enveloped and Rotavirus & Norovirus)	475	1	6		●
New 424105 ³	Tork Alcohol Gel Hand Sanitiser	S4	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Mycobacteria Virus	EN 1500, EN 12791, EN 13727, EN 1276, EN 1275, EN 13624, EN 1650, EN 14348, EN 14476 (all enveloped and Rotavirus & Norovirus)	1 000	1	6	●	●
511103	Tork Alcohol Gel Hand Sanitiser	Freestanding Bottle with Pump	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Mycobacteria Virus	EN 1500, EN 12791, EN 13727, EN 1276, EN 1275, EN 13624, EN 1650, EN 14348, EN 14476 (all enveloped and Rotavirus & Norovirus)	500	2	12		●
590103	Tork Alcohol Gel Hand Sanitiser	Pocket Bottle	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Mycobacteria Virus	EN 1500, EN 12791, EN 13727, EN 1276, EN 1275, EN 13624, EN 1650, EN 14348, EN 14476 (all enveloped and Rotavirus & Norovirus)	80		24		●
910103 ²	Tork Salubrin Hand Sanitiser	S1	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Virus	EN 1500, EN 1276, EN 1650, EN 14476 (all enveloped viruses)	1 000	1	6	●	
New 914103 ²	Tork Salubrin Hand Sanitiser	S4	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Virus	EN 1500, EN 1276, EN 1650, EN 14476 (all enveloped viruses)	1 000	1	6	●	
911103 ²	Tork Salubrin Hand Sanitiser	Freestanding Bottle with Pump	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Virus	EN 1500, EN 1276, EN 1650, EN 14476 (all enveloped viruses)	500	2	8		

¹ Use biocides safely. Always read the label and product information before use.

² Only available in UK.

³ Available from April 2021.

Tork Alcohol Hand Sanitisers

S1/S4, Freestanding Bottles, Accessories



420114



424114





511110



590110

Refills - Liquid¹

Art.No	Name	Refill ID	Quality	Fragrance	Colour	Efficacy	Efficacy Norms	Content/ Unit (ml)	Dose- Size (ml)	CON/ TRP		
420114	Tork Alcohol Liquid Hand Sanitiser	S1	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Mycobacteria Virus	EN 1500, EN 12791, EN 13727, EN 1276, EN 1275, EN 13624, EN 1650, EN 14348, EN 14476	1 000	1	6	●	●
New 424114 ²	Tork Alcohol Liquid Hand Sanitiser	S4	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Mycobacteria Virus	EN 1500, EN 12791, EN 13727, EN 1276, EN 1275, EN 13624, EN 1650, EN 14348, EN 14476	1 000	1	6	●	●
511110	Tork Alcohol Liquid Hand Sanitiser	Freestanding Bottle with Pump	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Mycobacteria Virus	EN 1500, EN 12791, EN 13727, EN 1276, EN 1275, EN 13624, EN 1650, EN 14348, EN 14476	500	2	12		●
590110	Tork Alcohol Liquid Hand Sanitiser	Pocket Bottle	Premium	Perfume Free	Dye Free	Bacteria Yeast & Fungi Mycobacteria Virus	EN 1500, EN 12791, EN 13727, EN 1276, EN 1275, EN 13624, EN 1650, EN 14348, EN 14476	80		24		●

Tork Alcohol-Free Hand Sanitisers



520202

Refills - Foam¹

Art.No	Name	Refill ID	Quality	Fragrance	Colour	Efficacy	Efficacy Norms	Content/ Unit (ml)	Dose-Size (ml)	CON/ TRP
520202-520402	Tork Alcohol-Free Foam Hand Sanitiser	S4	Premium	Perfume Free	Dye Free	Bacteria Yeast & Molds Virus	EN1500, EN1276, EN13727 EN EN1650, EN13624 EN 14476	1 000	1	6



511050



511051



511052



511058



511055



511054



511060



511080

Accessories

Art.No	Name	Refill ID	Colour	Material	Size HxWxD (mm)	CON/TRP
511050	Tork Bed Bracket for 500 ml Hand Sanitiser		Metal	Metal	110x79x114	20
511051	Tork Clip for 80 ml Hand Sanitiser		White	Plastic	75x35x5	50
511052	Tork Floor Protector for Skincare Dispenser	S1/S2/S4/S11	White	Metal	70x120x130	6
511058	Tork Floor Protector for Skincare Dispenser	S1/S2/S4/S11	Black	Metal	70x120x130	6
511055	Tork Hygiene Floor Stand	S1/S2/S4/F1	Anthracite	Metal	1 676x470x365	1
511054	AD-a-Glance® for Tork Hygiene Floor Stand		Black	Aluminium	400x125x125	1
511060	Tork Sanitiser Stand	S1/S2/S4	Black	Metal	1 321x500x500	1
511080	Tork Bracket for Armlever Dispenser	S1	White	Plastic	200x120	8

¹ Use biocides safely. Always read the label and product information before use.² Available from April 1st.