

Dispenser Key Chart

for Essity Professional Hygiene



No. 1 Key

Fits **current** Tork® dispensers:
70WM1, 83TR, and 87T

Fits obsolete Essity dispensers: 54, 54PW, 55, 55T, 55PW, 63C, 64, 64T, 65T, 66A, 67A, 70, 70A, 71, 71A, 71G, 73, 73P, 73PA, 73PW, 73T1, 74, 74C, 74CA, 74TA, 74T, 74W, 80WM, 81, 82, 82PW, 82T, 83B, 83W, 83T, 84, 84A, 84PW, 84T, 84W, 86T, 87A, 87PW, 92

Fits: Alwin, Bay West, Mosinee, Erving, Morgen, Timex, Stevens and Thompson, Nibroc, Marathon, Crown Zellerback, James River, Fort Howard



No. 200260 Key

Fits all **Tork Elevation** and **Tork Performance** dispensers

Tork Matic® Hand Towel Roll Dispenser System, Tork PeakServe® Continuous™ Hand Towel Dispenser System, Tork Soap and Sanitizer Dispenser Systems, Tork Jumbo, Twin Jumbo, Mini Jumbo, Twin Mini Jumbo, Coreless, SmartOne® and High Capacity Toilet Paper Roll Dispenser Systems, Tork Mini Centerfeed, Wiper Centerfeed, Tork Maxi Wiper Centerfeed, Tork Folded Wiper Cloth Dispenser and Tork Washstation Dispenser Systems.



No. 200261 Key

Key handle for all **Tork Elevation** and **Tork Performance** dispensers

Fits **obsolete** Tork Watertight Centerfeed dispenser:
302283



No. 470068 Key

Fits **Tork Coreless Bath Tissue** dispensers: 473000, 473008



No. 14 Key

Fits **current** Tork dispensers:
56TR, 59TR, 66TR, 67TR, 73TR,
84TR, 86ECO

Fits obsolete Tork dispensers: 52GR, 52T, 55TR, 56, 56A, 56PA, 56PW, 58, 58PW, 58T1, 65GR, 65T, 66PA, 66PW, 66T1, 68T, 75GR, 75T, 81T, 85GR, 85T, 86TR



Wausau Dual Key - 18002

Fits all Tork Electronic, Tork Mechanical, and Tork Mini-Mechanical dispensers; Tork High Capacity Dispensers for OptiCore®; and all Wausau Paper Dispensers.



No. 15 Key

Fits **current** Tork dispensers: 78T1,
93T, Sysco Complete 9209933

Fits **obsolete** Tork dispenser: 67TA

Fits Pre-1997 Wisconsin Tissue: 66