



Think ahead.

Tork Hand Sanitiser

**Premium quality hand sanitiser from the
#1 professional hygiene brand in the world**



**Highest rated
germ kill rate of
99.99%**





Don't compromise on hand hygiene - trust products from global leading brand Tork and feel confident with the security of ongoing local supply.

- **Maximum effectiveness** for sanitising hands against micro-organisms like bacteria and viruses with germ kill rate of 99.99%, complies with EN 13727
- **Premium formula** to protect hands with effective 80% alcohol and added moisturiser to maintain skin condition
- **User friendly** floor stand has optimum ergonomic dispenser height for all, including wheelchair users and good stability from 10.4kg steel base plate



420109
Tork Alcohol Gel
Hand Sanitiser
1,000mL



560000
Tork Skincare Dispenser
225 x 152 x 430mm



511056
Tork Tabletop Stand for
Skincare Dispensers
225 x 152 x 430mm



511057
Tork Tubular Floor
Hygiene Stand for
Skincare Dispensers
380 x 380 x 1370mm



Visit tork.com.au or tork.co.nz or contact your Tork Account Manager to order now



Think ahead.